

VISION V-4000 MICROSCOPE

HIGHLY RELIABLE AND COST-EFFECTIVE



The Vision V4000 binocular microscope is designed to fulfil all general microscopy requirements. It is a reliable, high-quality and very cost effective microscopy solution. The optical and mechanical parts of the microscope are manufactured to an exceptionally high standard, making this microscope high value for money.

Features:

- 30° inclined head rotates 360°
- Wide-field eyepieces - 10x/18mm
- Adjustable interpupillary distance on both eyepiece tubes
- Dioptic compensation on both eyepiece tubes
- DIN achromatic objectives (4x, 10x, 40x, 100x)
- Double layer mechanical stage, 125 x 130mm
- Coaxial coarse and fine focusing
- Cold-light white (fluorescent) 5W lamp
- N.A. Abbe 1.25 condenser with iris diaphragm
- Wide range voltage for international markets – 100 – 240V
- Plano-concave mirror
- **CE** certified
- Product code: WD6200

Precision.
Clarity.
Innovation.

Visit www.woodleylabdiagnostics.com for full product range