

## ▶ Refrigerated Centrifuge

available for hire or purchase for your clinical trial.



### Overview

Our refrigerated centrifuges are the ideal choice because of their lightweight and compact design making them ideal for international shipping.

They are easy to operate via a simple control panel, and the run can be programmed in either RPM or RCF - variable speed.

Various fixed angle and swing-out rotors are available to accommodate any kind of tube, the rotor to suit your needs can be supplied when placing your order for rental or purchase.

### Specification

- Rotor Options: A-8-17 rotor and A-4-38 rotor optional
- Max. Speed: 4400 RPM / 3000 x g
- Electronic Timer: 99 minutes
- Weight: 36 kg
- Dimensions: 260h x 380w x 580d mm
- Maintenance free induction motor
- Various capacity from 0.2ml to 85ml
- Voltage: 110V / 220V
- CE, UL and CAN/CSA certified

### Operating Features

- Adjustable: -9°C to +40°C
- Stand By Refrigeration: to maintain a particular temperature following use
- "Fast Temp" Function: to quickly reach the desired temperature
- "At Set RPM" Function: for precise centrifugation
- Two programming keys to store routine procedures
- Large dials and an easy-to-read display
- Switch between RPM/RCF
- Speed can be set from 100 RPM to a max. of 4,400 RPM, in increments of 100
- Separate "short-spin" button
- Timer can be set for up to 99 minutes
- Continuous mode can be selected
- Rotor options available:

Rotor	RPM	x g	Max. Capacity
A-8-17	4400	2800	8 x 15ml
A-4-38	4400	3000	4 x 85ml



*Please contact us for a quotation for the hire or purchase of this product.*

**Product description** - Eppendorf Refrigerated Centrifuge  
 Optional rotor A-817 or A-4-38

**Product code** - WD7CEN001/WD7CEN002



### UK Office

T: +44(0)8456 777001

F: +44(0)8456 777002

E: sales@woodleyequipment.com [www.woodleyequipment.com](http://www.woodleyequipment.com)

### USA Office

T: 1 800 471 9200 / 1 508 625 1693

F: 1 508 625 1721

**WOODLEY**  
 EQUIPMENT COMPANY LTD.