

# VISION V-4001 MICROSCOPE

## HIGHLY RELIABLE AND COST-EFFECTIVE



The Vision V4001 binocular microscope is designed to fulfil all general microscopy requirements. It is a reliable, high-quality and very cost effective microscopy solution. The optical and mechanical parts of the microscope are manufactured to an exceptionally high standard, making this microscope very high value for money. It comes complete with an internal rechargeable battery for use in the field.

### Features:

- 30° inclined head rotates 360°
- Wide-field eyepieces - 10x/18mm
- Adjustable interpupillary distance on both eyepiece tubes
- Variable light intensity
- Dioptic compensation on both eyepiece tubes
- DIN achromatic objectives (4x, 10x, 40x, 100x)
- Double layer mechanical stage, 125 x 130mm
- Coaxial coarse and fine focusing
- LED illuminator 50,000 hours use
- N.A. 1.25 Abbe condenser with iris diaphragm
- Wide range voltage for international markets – 100 – 240V
- Plano-concave mirror
- Battery pack: 3.6V NiMH 650 mA/h
- 8 hours battery charge time (microscope can be operated during charging)
- 7 hours battery life at maximum power
- CE certified
- Product code: WD6201



Visit [www.woodleylabdiagnostics.com](http://www.woodleylabdiagnostics.com) for full product range