

Clinical Thermometers

Glass clinical thermometer MERCURY FREE

Measuring range: +35°C +42 °C
Accuracy: +0,1°C / -0,15°C
Scale: paper
Prismatic capillary: blue
Non toxic, and environmentally friendly
Liquid : Ga, In, Sn
Plastic case.
Packed by 10 pieces.
€€1023



Order No. **801130**

EAN
8595131418395



Glass clinical thermometer MERCURY FREE with shaking case

Measuring range: +35°C +42 °C
Accuracy: +0,1°C / -0,15°C
Scale: paper
Prismatic capillary: blue
Non toxic, and environmentally friendly
Liquid : Ga, In, Sn
Plastic shaking case.
Packed in hang box.
€€1023



Order No. **801131**

EAN
8595131417695



Glass clinical thermometer with MERCURY

Measuring range: +35°C +42 °C
Accuracy: +0,1°C / -0,15°C
Scale: paper
Prismatic capillary: yellow
Liquid: Mercury
Plastic case
€€1014



Order No. **100001 C1**

EAN
8595131418081



Digital clinical thermometer EXATHERM Basic

Measuring range: +32,0°C +43,0°C
Accuracy : +/- 0,1 °C
Hard collar
Water Proof
Automatic switch off
Large display
Memory for last measurement
Fever alarm
€€0197



Order No. **700027**

EAN
8595131417671



Digital clinical thermometer EXATHERM

Measuring range: +32,0°C +43,0°C
Accuracy : +/- 0,1 °C
Elastic collar
Water Resistant
Automatic switch off
Large display
High Speed 30sec.
Memory for last measurement
Fever alarm
€€0197



Order No. **700023**

EAN
8595131418074



Digital clinical thermometer EXATHERM Rapid Comfort

Measuring range: +32,0°C +43,0°C
Accuracy : +/- 0,1 °C
Elastic collar
Water Resistant
Automatic switch off
Large display
High Speed 10sec.
Memory for last measurement
Fever alarm
€€0197



Order No. **700025**

EAN
8595131418401

